



NEW Thermo Scientific Sorvall RC 12BP Plus Large Capacity Centrifuge

**An integrated system for high-throughput,
large volume blood-banking and
bioprocessing applications**

The Best Just Keeps Getting Better.

The Thermo Scientific Sorvall RC 12BP Plus offers a complete processing system for high-throughput, large volume blood banking and bioprocessing applications.



An Integrated System Optimized for Your Success

- Sorvall® RC 12BP Plus centrifuge
- H-12000 Swinging Bucket rotor
- 2-liter True Fit Thermo Scientific Nalgene centrifuge bottles
- Array of adapters for process flexibility and convenience
- Thermo Scientific Centri-Log data collection software

An Integrated System Optimized for Your Success

Highest Throughput in a Single Step

- Process up to 12 blood bag systems in a single run – each with 500 mL
- Process up to 12 liters of volume per run

Versatile Accessory System

Supports an array of blood banking and bioprocessing applications and conveniently accommodates a broad range of blood bags, tubes, and bottles.

Simply Intuitive

- User-friendly interface makes set-up and operation intuitive and easy
- Enhanced display enables easy at-a-glance monitoring of set and run conditions
- 15-program memory provides one-touch recall of commonly used protocols

Ensure Run-to-Run Reproducibility

The Accumulated Centrifugal Effect™ (ACE) function assures reproducible blood product yields, regardless of rotor loading. ACE automatically compensates for variations in acceleration due to full or partial rotor loading by adjusting the centrifuge run time.

Maximize Blood Product Yields

In slow start/stop mode, select from 11 acceleration and deceleration profiles to customize results for maximum yields with minimal re-suspension.

Secure & Compliant

Centri-Log™ data collection software with security capabilities ensures total process control and run parameter traceability for cGMP, GMP, and GLP compliance.

Quiet & Comfortable

Quiet and compact with minimal heat output for comfortable, convenient in-laboratory placement. Low deck height facilitates loading and unloading of samples and cleaning of the unit. Support for standard electrical outlets provides added convenience.

Safety Engineered

A robust automated lid-lock control system protects users and makes lid openings and closures effortless.

2-Liter Nalgene® Polypropylene Bottles and Sealing Caps

Spin 2-liter Nalgene centrifuge bottles for highest capacities. Bottles are designed for an exact fit with the Thermo Scientific H-12000 rotor for optimal processing performance. A graduated wide-mouth simplifies decanting, re-suspensions, and pellet removal. Ergonomic design facilitates loading, recovery, and sample storage.

- Leakproof up to 7,340 xg
- Autoclavable
- Available in sterile form
- USP Class VI, non-cytotoxic and non-hemolytic



Thermo Scientific Sorvall centrifuges are the innovative choice for busy blood banks around the world.

Our large capacity centrifuges offer a unique set of benefits – dependability, maximum versatility, superior throughput and volume, along with simplicity, efficiency, and compliance.

Accelerate

your goals for performance and the bottom line.

- Reduce manpower, time, and energy costs through powerfully-efficient processing
- Scale up productivity through intuitive controls and paperless data collection/process document management
- Lower maintenance costs through well-proven Sorvall system dependability

Expand

your processing power.

- Process up to 12 x 450-500 mL blood bag systems simultaneously
- New adapters accommodate large numbers of blood tubes, allowing you to process up to 324 x 5/7 mL tubes at once
- A suite of proven accessories (swinging bucket rotor, 2 L bottles, and versatile adapter set) supports an array of blood banking and bioprocessing applications


Protect

your precious samples and your process development.

- Improve sample integrity through advanced capabilities – including Accumulated Centrifugal Effect (ACE) solutions
- Ensure cGMP compliance through easy-to-use Centri-Log data traceability software
- Reduce risk through well-proven Sorvall system dependability
- Security key restricts unauthorized access



State-of-the-art capabilities are accessed through a familiar, easy-to-use control panel. Frequently-used protocols can be recalled at the touch of a button to minimize time between runs.



**H-12000 Blood Processing Rotor,
Cat. No. 70050**

Proven Solutions, Optimized for Integrated Laboratory Workflows

The most reliable and versatile 12-liter low-speed centrifuge available. Compact and comfortable with operator-friendly controls, it is ideal for dependable, high-throughput sample processing of large volumes, at low speeds and constant temperatures.

- Maximum Speed: 4700 rpm
- Maximum RCF: 7340 xg
- Maximum Capacity: 12 L
- Spin 6 x 2000 mL Nalgene Bio-bottles
- Process 12 x 500 mL blood bag systems

Adapters and Carriers for Convenience & Versatility

A broad selection of adapters and carriers accommodate a variety of blood bags, tubes, and bottles. Large capacity tube accommodation with select adapters.

Centri-Log Data Collection Software

Electronic data capture software for the documentation (recording and monitoring) of centrifugation processes.

Step-Run Mode Automation

Buffy coat protocols can be optimized using the step-run mode, which automatically increments run speed in pre-programmed stages. Customize results for maximum yields with minimum suspension and maximized blood product yields.

Advanced Capabilities for Demanding Applications

Centri-Log Data Collection Software

A paperless process tracking and compliance solution for large capacity, low-speed centrifuge systems.

This innovative software provides effective documentation of centrifugation processes, enabling lifecycle management of processed samples, equipment optimization, and compliance with standard operating procedures.

Centri-Log Viewer:

Program File Evaluation

Centrifuge type
RC12BP+

User loading ☐

User unloading ☐

Product ☐

Blood bags

408732842
408732849
408732839
408732844

Run status ☐

No monitoring error



20100623_100959_01.CCP

Help

Centrifuge identification

Production #1

Time of process

23.06.2010

10:09:00

7842

Dr. Tansey

7842

Dr. Tansey

WHBD

Blood processing step #1


Process defined for

RC12BP+

Definition / set values

Definition	Set values
8	8
4	4
3620	3620
4354	4354
12:00	12:00
20	
H12000	

Graphics




Print



Close

Tolerance limits



Print




Close

- Quick, easy, secure, and accurate data collection and document management
- Continual monitoring of speed, time, and temperature with alarm messages in case of deviation from protocol
- Ensured sample security and integrity
- Ability to check run parameters against procedures in PC database
- Connect to central database through customized export file
- Eliminate manual recording errors
- Reduce risk of equipment downtime and quality failures
- Centri-Log data collection software network can be customized to fit the needs of your laboratory
- Support for standard technologies: Windows® compatible with RS-232 port/USB: Windows XP and earlier versions: (Windows NT, Windows 98, Windows ME, and Windows 2000)

The Sorvall RC 12BP Plus Centrifuge System is the ideal solution for large-batch bioprocessing applications in biotechnology and pharmaceutical laboratories.





H-12000 Bioprocessing Rotor,
Cat. No. 70080

Powering the most demanding processing applications

Process Bacterial and Yeast Cultures and Clarify Broths


A 6 x 2000 mL bottle capacity allows you to process a 12-liter bioreactor harvest in a single run. The ACE function and slow start/stop modes maximize yields and reproducibility when processing multiple batches.

Harvest Mammalian Cell and Tissue Cultures

Adapters for the H-12000 rotor simplify cell and tissue culture harvests, letting you spin sample vessels ranging from 2 L bottles down to 50/15 mL conical tubes.

Protect Valuable Samples

With the Sorvall RC 12BP Plus Centrifuge System, you can control acceleration and deceleration rates for delicate low-speed pelleting applications.



Thermo Scientific Rotors and Adapters

Cat. No.	Description	Rotor Capacity (Places x Volume)	Max Speed (rpm)	Max RCF (xg)
70050	H-12000 Blood Processing Swinging Bucket Rotor (includes 6 buckets)	12 x 500 mL		
	Adapters and Accessories for Blood Processing Rotor			
77030	① 450-500 mL quads with additive	12 x 500 mL	4,700	7,340
77032	② 450-500 mL quads without additive	12 x 500 mL	4,700	7,340
77031	③ 200-250 mL platelet rich plasma or buffy coat	12 x 300 mL	4,700	7,340
77033	Covers for adapters			
77065	Adapter rack (6 place)			
77064	Adapter rack (single)			
75007065	④ BIOLink Adapter for H-12000		4,000	4,613
75007058	⑤ 5/7 mL tube*	6 x 49 = 294 mL	4,000	4,613
75007057	⑥ 10 mL tube*	6 x 37 = 222 mL	4,000	4,613
75007056	⑦ 15 mL conical tube*	6 x 22 = 132 mL	4,000	4,613
75007055	⑧ 50 mL conical tube*	6 x 10 = 60 mL	4,000	4,613
ME060001	⑨ 5/7 mL blood collection tube (Greiner®)	6 x 54 = 324 mL	4,000	4,613
ME080001	⑩ 5/7 mL blood collection tube (BD®)	6 x 54 = 324 mL	4,000	4,613
*Note: Adapters in lines 5 through 8 require the BIOLink Adapter.				
70080	H-12000 Bioprocessing Swinging Bucket Rotor (includes 6 buckets)	6 x 2000 mL		
	Adapters for Bioprocessing Rotor			
3120-2002	2 Liter Nalgene Bio-bottle	6 x 2000 mL	4,700	7,340
3120-1000	1 Liter polypropylene bottle	6 x 1000 mL	4,700	7,340
3122-1000	1 Liter polycarbonate bottle	6 x 1000 mL	4,700	7,340
77074 1	Liter bottle adapter	6 x 1000 mL	4,700	7,340

New adapter capabilities for the H-12000 rotor maximize process performance and application versatility.

Conical Tube

- 10 x 50 mL conical or 60 per run

Blood Collection Tube

- 54 x 5/7 mL blood collection tubes or 324 per run



Specifications and Ordering Information

Sorvall RC 12BP Plus Centrifuge System

Centrifuge Specifications

Maximum speed	4,700 rpm
Maximum RCF	7,340 xg
Maximum capacity	12 L, 6 x 2000 mL Nalgene Bio-bottles, 12-500 mL quad blood bags
Drive	High-performance induction
Imbalance detector	50 gram tolerance
Accel/decel profiles	11
Run time	1 min. - 99 min 59 sec., continuous operation
Temperature control range	-10 to +40 °C
Dimensions (H x W x D)	1178 x 800 x 905 mm (46.3 x 31.5 x 35.6 in)
Weight	470 kg (1,036 lbs)
Control	Microprocessor
Programmability	15 programs
Functions	Centri-Log traceability software, ACE, Step-mode runs, slow start/stop
Safety features	Automatic lid lock, imbalance detection, armored guard ring, security key switch
Certification	cULus approved and CE marked

Ordering Information

Description	Cat. No.
Sorvall RC 12BP Plus, Single Phase, 208-220V, 60Hz	75007030
Sorvall RC 12BP Plus, Single Phase, 230V, 60Hz	75007031
Sorvall RC 12BP Plus, Single Phase, 220V-230V, 50Hz	75007032
Sorvall RC 12BP Plus, 3Phase, 400V, 50Hz	75007033

Centri-Log Data Collection Software

Windows compatible with RS-232 port/USB: Windows XP and earlier versions: (Windows NT, Windows 98, Windows ME, and Windows 2000)
Computer not included with software. Customer to provide PC.

Description	Cat. No.
Basic kit for PC w/USB (required for system setup)	75007535
Networking kit, includes bar code reader (kit required for each centrifuge networked)	75007802
Installation kit for first centrifuge (includes power supply)	75007516
Installation kit for each additional centrifuge	75007517

Contact your sales representative to discuss the right Centri-Log networking solution for your laboratory.

Comprehensive Centrifuge Solutions

We offer a wide range of Sorvall centrifuges for the pursuit of scientific discoveries. As the world leader in serving science, we also offer a full range of high-quality solutions for your laboratory. Visit www.thermoscientific.com for more information. To learn about our full range of Thermo Scientific centrifuge solutions, including benchtop models and microcentrifuges, visit www.thermoscientific.com/centrifuge.

Worldwide Service and Support

We are committed to keeping your lab equipment working at peak performance levels. Our goal is to help you lower ownership costs, manage labs more effectively, and increase productivity. Contact your Thermo Scientific product representative to learn more about our service offerings, including service agreements, preventative maintenance, onsite field repair, depot repair, compliance services, and educational services.



© 2010 Thermo Fisher Scientific Inc. All rights reserved. Windows is a registered trademark of Microsoft Corporation in the United States and other countries. BD is a registered trademark of Beckton, Dickinson and Company. Greiner is a trademark of Greiner Bio-One International. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Validation Services for Peace of Mind

Our certified validation engineers will ensure that your Sorvall RC 12BP Plus Centrifuge System meets designated compliance criteria at the time of installation and throughout its long operational life. Our validation services also include documentation support to address today's stringent regulatory requirements. We offer a suite of packages that ensure compliance with applicable Installation, Operational and Performance Qualification standards (IQ/OQ/PQ), and can develop customized plans to meet your specific laboratory needs.

North America: USA/Canada +1 866 984 3766 (866-9-THERMO)

www.thermoscientific.com/centrifuge

Europe: Austria +43 1 801 40 0, Belgium +32 2 482 30 30, France +33 2 2803 2180, Germany national toll free 08001-536 376, Germany international +49 6184 90 6940, Italy +39 02 02 95059 434-254-375, Netherlands +31 76 571 4440, Nordic countries +358 9 329 100, Russia/CIS +7 (812) 703 42 15, Spain/Portugal +34 93 223 09 18, Switzerland +41 44 454 12 12, UK/Ireland +44 870 609 9203

Asia: China +86 21 6865 4588 or +86 10 8419 3588, India toll free 1800 22 8374, India +91 22 6716 2200, Japan +81 45 453 9220,

Other Asian countries +852 2885 4613, **Countries not listed:** +49 6184 90 6940 or +33 2 2803 2000

BRCFGRC12BPPLUS 1110

Thermo
SCIENTIFIC