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REFERÊNCIAS

1. Rudolph J, Marcantonio E. Postoperative delirium: Acute change with long-term implications. *Anesth Anal.* 2011;112(5):1202–1211.
2. Wong J, Song D, Blanshard H, et al. Titration of isoflurane using BIS index improves early recovery of elderly patients undergoing orthopedic surgeries. *Can J Anaesth.* 2002;49(1):13–18.
3. Sebel PS, Bowdle TA, Ghoneim MM, et al. The incidence of awareness during anesthesia: a multicenter United States study. *Anesth Analg.* 2004;99(3):833–839.
4. Luginbuhl M, Wuthrich S, Petersen-Felix S, Binden A, Schnider T. Different benefit of bispectral index (BIS™) in desflurane and propofol anesthesia. *Acta Anaesthesiol Scand.* 2003;47(2):165–173.
5. Gan T, Glass P, Windsor A, et al. Bispectral index monitoring allows faster emergence and improved recovery from propofol, alfentanil, and nitrous oxide anesthesia. *Anesthesiology.* 1997;87(4):808–815.
6. Punjasawadwong Y, Phongchiewboon A, Bunchungmongkol N. Bispectral index for improving anaesthetic delivery and postoperative recovery (Review). *Cochrane Database Syst. Rev.* 2014;(6):1–11.
7. Chan M, Cheng B, Lee T, Gin T, CODA Trial Group. BIS-guided anesthesia decreases postoperative delirium and cognitive decline. *J Neurosurg Anesthesiol.* 2013;25(1):33–42.
8. Lewis SR, Pritchard MW, Fawcett LJ, Punjasawadwong Y. Bispectral index for improving intraoperative awareness and early postoperative recovery in adults. *Cochrane Database Syst. Rev.* 2019;26(9):CD003843.
9. Punjasawadwong Y, Chau-In W, Laopaiboon M, Punjasawadwong S, Pin-On P. Processed electroencephalogram and evoked potential techniques for amelioration of postoperative delirium and cognitive dysfunction following non-cardiac and non-neurosurgical procedures in adults. *Cochrane Database Syst. Rev.* 2018;5:CD01128.
10. Liu SS. Effects of Bispectral index monitoring on ambulatory anesthesia: a meta-analysis of randomized controlled trials and cost analysis. *Anesthesiology.* 2004;101(2):311–5.
11. Song D, Joshi GP, White PF. Titration of volatile anesthetics using bispectral index facilitates recovery after ambulatory anesthesia. *Anesthesiology.* 1997;87(4):842–8.
12. Myles P, Leslie K, McNeil J, Forbes A, Chan M. Bispectral index monitoring to prevent awareness during anaesthesia: the B-Aware randomised controlled trial. *Lancet.* 2004;363(9423):1757–1763.
13. Ekman A, Lindholm M, Lennmarken C, Sandin R. Reduction in the incidence of awareness using BIS monitoring. *Acta Anaesthesiol Scand.* 2004;48(1):20–26.
14. Radtke F, Franck M, Lendner J, Kruger S, Wernecke K, Spies C. Monitoring depth of anaesthesia in a randomized trial decreases the rate of postoperative delirium but not postoperative cognitive dysfunction. *British Journal of Anaesthesia.* 2013;110(S1):98–105.
15. Sieber F, Zakriya K, Gottschalk A, et al. Sedation depth during spinal anesthesia and the development of postoperative delirium in elderly patients undergoing hip fracture repair. *Mayo Clin Proc.* 2010;85(1):18–26.
16. Klopman M, Sebel P. Cost-effectiveness of bispectral index monitoring. *Current Opinion in Anesthesiology.* 2011;24(1):177–181.
17. Fritz B, Kalarickal P, Maybrier H, et al. Intraoperative electroencephalogram suppression predicts postoperative delirium. *International Anesthesia Research Society.* 2015;1–9. doi: 10.1213/ANE.0000000000000989.
18. Whitlock E, Torres B, Lin N, et al. Postoperative delirium in a substudy of cardiothoracic surgical patients in the bag-recall clinical trial. *Neuroscience in Anesthesiology and Perioperative Medicine.* 2014;118(4):809–817.
19. Liu J, Singh H, White P. Electroencephalographic bispectral index correlates with intraoperative recall and depth of propofol-induced sedation. *Anesth Analg.* 1997;84(1):185–189.

REFERÊNCIAS

20. Ahmad S, Yilmaz M, Marcus RJ, Glisson S, Kinsella A. Impact of bispectral index monitoring on fast tracking of gynecologic patients undergoing laparoscopic surgery. *Anesthesiology*. 2003 Apr; 98(4): 849–852.
21. Satisha M, Sanders GM, Badrinath MR, Ringer JM, Morley AP. Introduction of bispectral index monitoring in a district general hospital operating suite: a prospective audit of clinical and economic effects. *Eur J Anaesthesiol*. 2010 27(2): 196–201.
22. Shepherd J, Jones J, Frampton G, Bryant J, Baxter L, Cooper K. Clinical effectiveness and cost effectiveness of depth of anaesthesia monitoring (E-Entropy, Bispectral Index and Narcotrend): a systematic review and economic evaluation. *Health technology assessment (Winchester, England)*. 2013; 17(34): 1–264.
23. Zhang C, Xu L, Ma Y-Q, et al. Bispectral index monitoring prevent awareness during total intravenous anesthesia: a prospective, randomized, double-blinded, multi-center controlled trial. *Chin Med J (Engl)*. 2011; 124(22): 3664–9.
24. Liu S. Effects of bispectral index monitoring on ambulatory anesthesia. *Anesthesiology*. 2004; 101(2): 311–315.
25. Glass P, Bloom M, Kears L, Rosow C, Sebel P, Manberg P. Bispectral analysis measures sedation and memory effects of propofol, midazolam, isoflurane, and alfentanil in healthy volunteers. *Anesthesiology*. 1997; 86(4): 836–847.
26. Brown D, Avramov M, Tuman K, Ivankovich A. The effect of dexmedetomidine on EEG Bispectral index. *Anesthesia and Analgesia*. 1999; 88(2S): 52S.
27. Turkmen A, Altan A, Turgut N, Vatansever S, Gokkaya S. The correlation between the Richmond agitation-sedation scale and bispectral index during dexmedetomidine sedation. *European Journal of Anaesthesiology*. 2006; 23(4): 300–304.
28. Flash R, Windsor A, Sigl J, Sebel P. Recovery of consciousness after thiopental or propofol. *Anesthesiology*. 1997; 86(3): 613–619.
29. Leslie K, Sessler D, Schroeder M, Walters K. Propofol blood concentration and the bispectral index predicts suppression of learning during propofol/epidural anesthesia in volunteers. *Anesth Analg*. 1995; 81(6): 1269–1274.
30. Tripathi M, Kumar V, Kalashetty M, Malviya D, Bais P, Sanjeev O. Comparison of dexmedetomidine and midazolam for sedation in mechanically ventilated patients guided by bispectral index and sedation-agitation scale. *Anesth Essays Res*. 2017; 11(4): 828–833.
31. Dahaba A. Different condition that could result in the bispectral index indicating an incorrect hypnotic state. *Anesth Analg*. 2005; 101(3): 765–773.
32. BIS™ 1.0 Owners Manual. PN: PT00087175X01. 2018.
33. Denman WT, Swanson EL, Rosow D, Ezbicki K, Connors PD, Rosow CE. Pediatric evaluation of the Bispectral index (BIS™) monitor and correlation of BIS™ with end-tidal sevoflurane concentration in infants and children. *Anesth Analg*. 2000; 90(4): 872–7.
34. Davidson AJ, McCann ME, Devavaram P, Auble SA, Sullivan LJ, Gillis JM, Laussen PC. The differences in the Bispectral index between infants and children during emergence from anesthesia after circumcision surgery. *Anesth Analg*. 2001; 93(2): 326–30.
35. McCann ME, Bacsik J, Davidson A, Auble S, Sullivan L, Laussen P. The correlation of Bispectral index with end-tidal sevoflurane concentration and haemodynamic parameter in preschoolers. *Paediatr Anaesth*. 2002; 12(6): 519–25.
36. Whyte SD, Booker PD. Bispectral index during isoflurane anesthesia in pediatric patients. *Anesth Analg*. 2004; 98(6): 1644–9.
37. Degoute CS, Macabeo C, Dubreuil C, Duclaux R, Bannillon V. EEG Bispectral index and hypnotic component of anaesthesia induced by sevoflurane: comparison between children and adults. *Br J Anaesth*. 2001; 86(2): 209–12.
38. Bannister CF, Brosius KK, Sigl JC, Meyer BJ, Sebel PS. The effect of Bispectral index monitoring on anesthetic use and recovery in children anesthetized with sevoflurane in nitrous oxide. *Anesth Analg*. 2001; 92(4): 877–81.

REFERÊNCIAS

39. Purdon P, Sampson A, Pavone K, Brown E. Clinical electroencephalography for anesthesiologists. *Anesthesiology*. 2015;123(4):937–960.
40. Jameson L, Sloan T. Using EEG to monitor anesthesia drug effects during surgery. *Journal of Clinical Monitoring and Computing*. 2006;20:445–472.
41. Choi B. Characteristics of electroencephalogram signatures in sedated patients induced by various anesthetic agents. *J Dent Anesth Pain Med*. 2017;17(4):241–251.
42. George E, Devadas G, Avinash K, Vijayan S. Study of different spectral parameters to measure depth of anesthesia. *Biomedical Research*. 2017;28(15):6533–6540.
43. Mashour GA, Avidan MS. Black swans: challenging the relationship of anaesthetic-induced unconsciousness and electroencephalographic oscillations in the frontal cortex. *Brit J Anaesthesia*. 2017;119(4):563–565.
44. Hagihara S. Changes in the electroencephalogram during anaesthesia and their physiological basis. *British Journal of Anaesthesia*. 2015;27–31.
45. Hotz MA, Ritz R, Linder L. Auditory and electroencephalographic effects of midazolam and a hydroxy-midazolam in healthy subjects. *Br J Clin Pharmacol*. Jan 2009;49(1):72–79.
46. Hajat Z, Ahmad N, Andrzejewski J. The role and limitations of EEG-based depth of anaesthesia monitoring in theatres and intensive care. *Anaesthesia*. 2017;72(1):38–47.
47. Purdon P, Pavone K, Akeju O, et al. The Ageing Brain: Age-dependent changes in the electroencephalogram during propofol and sevoflurane general anaesthesia. *British Journal of Anaesthesia*. 2015;46–57.

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