

# [ Outstanding Performance ]

Across a Broader Range of Tissues and Applications

## Endo GIA™ Reloads with Tri-Staple™ Technology

Covidien's Tri-Staple™ Technology provides Endo GIA™ Reloads with outstanding performance across a broader range of tissue thicknesses, thereby simplifying the cartridge selection process.



### Outstanding Performance

- Improved burst pressure strength\*
- Reduced tissue compression stress on outer staple row compared to Endo GIA™ Universal reloads
- Consistent hemostatic performance across a broader indicated tissue range
- Improved tissue retention during manipulation and transection

### Broader Range of Tissues and Applications

- Broader indicated tissue thickness ranges for tan, purple and black reloads
- Unique black reloads for extra-thick tissue applications previously beyond the indications of any MIS stapler
- Curved tip reload enhances visualization, blunt dissection and tissue manipulation
- Medium/thick tissue purple reload can be introduced through a 12mm port

### Simplification of Cartridge Selection

RELOAD SELECTION CHART				
Extra-Thin Range	Thin Range	Medium Range	Thick Range	Extra-Thick Range
Vascular				
Endo GIA™ Reloads without Tri-Staple™ Technology				
2.0mm	2.5mm	3.5mm	4.8mm	
12mm port		15mm port		
New Endo GIA™ Reloads with Tri-Staple™ Technology				
Extra-Thin/Vascular*				
		Vascular/Medium		
		Medium/Thick		Extra-Thick
		12mm port	15mm port	
Ethicon Echelon™/ETS Reloads				
Thin/Mesenteric	Vascular/Thin	Regular	Regular/Thick	Thick
		12mm port		

\* The results from bench testing may not directly correlate with clinical performance.

† Extra-Thin/Vascular reload contains three rows of 2.0mm staples (non-varied-height) and contains all of the other features and benefits of an Endo GIA™ Reload with Tri-Staple™ Technology.


TM\* Trademark of its respective owner.




**COVIDIEN**

*positive results for life™*

## Endo GIA™ Reloads with Tri-Staple™ Technology



### Varied-height Staples and Stepped Cartridge Face



\*Products depicted are not actual size.

### Ordering Information

Reorder Code	Description	Cartridge Color	Staple Size (inner to outer row)
EGIA45AV*	Endo GIA™ 45mm Articulating Vascular Reload	Gray	2mm, 2mm, 2mm
EGIA45CTAV*	Endo GIA™ 45mm Curved Tip Articulating Vascular Reload	Gray	2mm, 2mm, 2mm
EGIA30AVM	Endo GIA™ 30mm Articulating Vascular/Medium Reload with Tri-Staple™ Technology	Tan	2mm, 2.5mm, 3mm
EGIA30CTAVM	Endo GIA™ 30mm Curved Tip Articulating Vascular/Medium Reload with Tri-Staple™ Technology	Tan	2mm, 2.5mm, 3mm
EGIA45AVM	Endo GIA™ 45mm Articulating Vascular/Medium Reload with Tri-Staple™ Technology	Tan	2mm, 2.5mm, 3mm
EGIA45CTAVM	Endo GIA™ 45mm Curved Tip Articulating Vascular/Medium Reload with Tri-Staple™ Technology	Tan	2mm, 2.5mm, 3mm
EGIA60AVM	Endo GIA™ 60mm Articulating Vascular/Medium Reload with Tri-Staple™ Technology	Tan	2mm, 2.5mm, 3mm
EGIA60CTAVM	Endo GIA™ 60mm Curved Tip Articulating Vascular/Medium Reload with Tri-Staple™ Technology	Tan	2mm, 2.5mm, 3mm
EGIA30AMT	Endo GIA™ 30mm Articulating Medium/Thick Reload with Tri-Staple™ Technology	Purple	3mm, 3.5mm, 4mm
EGIA45AMT	Endo GIA™ 45mm Articulating Medium/Thick Reload with Tri-Staple™ Technology	Purple	3mm, 3.5mm, 4mm
EGIA60AMT	Endo GIA™ 60mm Articulating Medium/Thick Reload with Tri-Staple™ Technology	Purple	3mm, 3.5mm, 4mm
EGIA60CTAMT	Endo GIA™ 60mm Curved Tip Articulating Medium/Thick Reload with Tri-Staple™ Technology	Purple	3mm, 3.5mm, 4mm
EGIA45AXT	Endo GIA™ 45mm Articulating Extra-Thick Reload with Tri-Staple™ Technology	Black	4mm, 4.5mm, 5mm
EGIA60AXT	Endo GIA™ 60mm Articulating Extra-Thick Reload with Tri-Staple™ Technology	Black	4mm, 4.5mm, 5mm

\*Reorder codes EGIA45AV and EGIA45CTAV, while offering a slimmer fixed anvil, are vascular reloads that do not come with Tri-Staple™ Technology.

### IMPORTANT

Please refer to the package insert for complete instructions, contraindications, warnings and precautions.



COVIDIEN, COVIDIEN with logo, Covidien logo and *positive results for life* are U.S. and/or internationally registered trademarks of Covidien AG.  
© 2010 Covidien.

06/10 P100159

For more information or to set up a product demonstration, contact your local Covidien Surgical Device Representative or call Customer Service at 1-800-722-8772.

60 MIDDLETOWN AVENUE  
NORTH HAVEN, CT  
06473

1-800-722-8772 [T]  
1-888-636-1002 [F]

WWW.TRISTAPLE.COM