



QIAstat-Dx Respiratory Panel V2



The next generation of syndromic insights

The QIAstat-Dx™ (DiagCORE®) Respiratory Panel V2 utilizes proven QIAGEN® sample and assay technologies to provide clinical insights for patients with suspected upper respiratory infections.



Bacterial

- *Mycoplasma pneumoniae*
- *Legionella pneumophila*
- *Bordetella pertussis*

Viral

- Influenza A
- Influenza A subtype H1N1/2009
- Influenza A subtype H1
- Influenza A subtype H3
- Influenza B
- Coronavirus 229E
- Coronavirus HKU1
- Coronavirus NL63
- Coronavirus OC43
- Parainfluenza virus 1
- Parainfluenza virus 2
- Parainfluenza virus 3
- Parainfluenza virus 4
- Respiratory Syncytial virus A/B
- Human Metapneumovirus A/B
- Adenovirus
- Bocavirus
- Rhinovirus/Enterovirus

Ordering Information - Cartridge

Product	Contents	Product code
QIAstat-Dx (DiagCORE) Respiratory Panel V2	6 individually packaged cartridges containing all reagents needed for sample preparation and multiplex RT-real time PCR plus internal control, including 6 transfer pipettes	210002

Ordering Information - Instrument

Product	Contents	Product code
QIAstat-Dx (DiagCORE) Analyzer 1.0	Instrument consists of 1 each of both the Operational Module and Analytical Module.	100000
QIAstat-Dx (DiagCORE) Analytical Module	1 each of module containing hardware and software for sample testing and analysis	110001
QIAstat-Dx (DiagCORE) Operational Module	1 each of module to enable interaction with the Analyzer	110002

The QIAstat-Dx (DiagCORE) is intended for in vitro diagnostic use in Europe. Not available in all countries.

Discover QIAstat-Dx – the next generation of syndromic insights

Visit **QIAstat-Dx.com** for more info

Trademarks: QIAGEN®, Sample to Insight®, QIAstat-DX™ (QIAGEN Group). DiagCORE® (STAT-Dx Life S.L.) Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

1114003 06/2018 PROM-12460-001 © QIAGEN 2018, all rights reserved.

Ordering www.qiagen.com/shop | Technical Support support.qiagen.com | Website www.qiagen.com