

## cell culture promotions

improving life in your lab



# grow spheroid cultures with reproducible results

## Thermo Scientific<sup>™</sup> Nunclon<sup>™</sup> Sphera Flasks, **Dishes and Plates**

Buy 2 cases of Nunclon Sphera treated flasks, dishes or plates and receive a third case at no additional charge.



Qualifying Products for Promotion			
Cat. No.	Description	Qty/Bag	Qty/Case
174925	Microwell 96U-Well Plate, Round Bottom, 300µL, sterile; with lid	1	8
174927	Microwell 96U-Well Plate, Flat Bottom, 400µL, sterile; with lid	1	8
174930	Multidish 24-Well, Culture Area 1.9cm <sup>2</sup> , suggested working volume 1mL, sterile; with lid	1	7
174931	Multidish 12-Well, Culture Area 3.8cm <sup>2</sup> , suggested working volume 2mL, sterile; with lid	1	7
174932	Multidish 6-Well, Culture Area 9.6cm <sup>2</sup> , suggested working volume 3mL, sterile; with lid	1	7
174943	Dish 35 MM, Culture Area 8.8cm <sup>2</sup> , suggested working volume 3mL, sterile; with lid	5	20
174944	Dish 60 MM, Culture Area 221.5cm <sup>2</sup> , suggested working volume 5mL, sterile; with lid	5	20
174945	Dish 90 MM, Culture Area 256.7 cm <sup>2</sup> , suggested working volume 12.5 mL, sterile; with lid	5	20
174951	T25 Cell Culture Flask, Culture Area 25cm <sup>2</sup> , suggested working volume 7mL, sterile; with cap	6	18
174952	T25 Cell Culture Flask, Culture Area 75cm <sup>2</sup> , suggested working volume 7mL, sterile; with cap	4	24

Learn more about Nunclon Sphera and find supporting.

#### last line of defense

against contamination

## Thermo Scientific<sup>™</sup> Nalgene Rapid-Flow<sup>™</sup> Filters

Buy 2 cases of Rapid-Flow filters and get a third case at no additional charge.



Qualifying Products for Promotion		
Cat. No.	Description	Qty/Case
565-0010	Filter Unit, 150mL, 0.1µm aPES Membrane, 50mm diameter	12
568-0010	Filter Unit, 250mL, 0.1µm aPES Membrane, 50mm diameter	12
566-0010	Filter Unit, 500mL, 0.1µm aPES Membrane, 75mm diameter	12
567-0010	Filter Unit, 1000mL, 0.1µm aPES Membrane, 90mm diameter	12
564-0020	Filter Unit, 50mL, 0.2µm aPES Membrane, 50mm diameter	12
565-0020	Filter Unit, 150mL, 0.2µm aPES Membrane, 50mm diameter	12
568-0020	Filter Unit, 250mL, 0.2µm aPES Membrane, 75mm diameter	12
566-0020	Filter Unit, 500mL, 0.2µm aPES Membrane, 90mm diameter	12
569-0020	Filter Unit, 1000mL, 0.2µm aPES Membrane, 90mm diameter	12
567-0020	Filter Unit, 1000mL, 0.2µm aPES Membrane, 90mm diameter	12
596-3320	Bottle Top Filter, 150mL, 0.2µm aPES Membrane, 50mm diameter, 33mm Neck	12
596-4520	Bottle Top Filter, 150mL, 0.2µm aPES Membrane, 50mm diameter, 45mm Neck	12
595-3320	Bottle Top Filter, 500mL, 0.2µm aPES Membrane, 50mm diameter, 33mm Neck	12
595-4520	Bottle Top Filter, 500mL, 0.2µm aPES Membrane, 75mm diameter, 45mm Neck	12
597-3320	Bottle Top Filter, 1000mL, 0.2µm aPES Membrane, 90mm diameter, 33mm Neck	12
597-4520	Bottle Top Filter, 1000mL, 0.2µm aPES Membrane, 90mm diameter, 45mm Neck	12

Learn more about our Rapid-Flow filters for cell culture.

### small filters

big results



## Thermo Scientific<sup>™</sup> Nalgene Syringe Filters

Buy 2 cases of sterile syringe filters and get a third case at no additional charge.

Qualifying Products for Promotion		
Cat. No.	Description	Qty/Case
720-1320	13mm diameter, 2 to 10mL sample size, PES Membrane, 0.2µm pore size, sterile	50
725-2520	25mm diameter, 10 to 50mL sample size, PES Membrane, 0.2 $\mu$ m pore size, sterile	50
725-2545	25mm diameter, 10 to 50mL sample size, PES Membrane, 0.45 $\mu$ m pore size, sterile	50
723-2520	25mm diameter, 10 to 50mL sample size, SFCA Membrane, $0.2\mu m$ pore size, sterile	50
723-2545	25mm diameter, 10 to 50mL sample size, SFCA Membrane, $0.45 \mu m$ pore size, sterile	50
726-2520	25mm diameter, 10 to 50mL sample size, Nylon Membrane, $0.2\mu m$ pore size, sterile	50
726-2545	25mm diameter, 10 to 50mL sample size, Nylon Membrane, 0.45µm pore size, sterile	50

Learn more about our syringe filter products.

## innovation

at every turn



## Thermo Scientific<sup>™</sup> Nunc<sup>™</sup> Conical Tubes

Buy 2 cases of 15mL, 50mL or 250mL conical tubes and a third case at no additional charge.

Qualifying Products for Promotion		
Cat. No.	Description	Qty/Case
339650	Conical Tubes, 15mL, Bulk Packaged	50/500
339651	Conical Tubes, 15mL, Racked	25/500
339652	Conical Tubes, 50mL, Bulk Packaged	25/500
339653	Conical Tubes, 50mL, Bulk Packaged, Racked	25/300
376814	Conical Tubes, 250mL, Bagged	4/40

Learn more about our conical tubes.

#### superior precision

and performance

## Thermo Scientific<sup>™</sup> Nunc<sup>™</sup> Serological Pipettes

Buy 2 cases of serological pipettes and receive a third case at no additional charge.



#### **Qualifying Products for Promotion** Cat. No. Serological Pipettes, Individually wrapped, Paper Peel Packaging 170353 1mL Pipette, Individually wrapped, paper peel packaging, plugged 100/1000 170354 2mL Pipette, Individually wrapped, paper peel packaging, plugged 100/500 170355 50/200 5mL Pipette, Individually wrapped, paper peel packaging, plugged 170356 10mL Pipette, Individually wrapped, paper peel packaging, plugged 50/200 170357 25mL Pipette, Individually wrapped, paper peel packaging, plugged 50/200 170358 25/100 50mL Pipette, Individually wrapped, paper peel packaging, plugged Shortie, Individually wrapped, Paper Peel Packaging 170360 5mL shortie pipette 50/200 170361 10mL shortie pipette 50/200 Serological Pipettes, Individually wrapped, Plastic Film Packaging 170364 1mL Pipette, Individually wrapped, plastic film packaging, plugged 100/1000 170365 2mL Pipette, Individually wrapped, plastic film packaging, plugged 100/500 170366 5mL Pipette, Individually wrapped, plastic film packaging, plugged 50/200 170367 10mL Pipette, Individually wrapped, plastic film packaging, plugged 50/200 170368 50/200 25mL Pipette, Individually wrapped, plastic film packaging, plugged 170369 50mL Pipette, Individually wrapped, plastic film packaging, plugged 25/100 Serological Pipettes, Bulk Packed 170371 100/1000 1mL Pipette, bulk packed, plugged 170372 2mL Pipette, bulk packed, plugged 100/500 170373 50/200 5mL Pipette, bulk packed, plugged 170374 50/200 10mL Pipette, bulk packed, plugged 170375 50/200 25mL Pipette, bulk packed, plugged 170376 50mL Pipette, bulk packed, plugged 25/100

Learn more about our serological pipettes.

# engineered for superior and consistent cell growth

## **Thermo Scientific<sup>™</sup>Nunc<sup>™</sup>** Flasks, Dishes & Multidishes

Buy 2 cases of Nunclon™ Delta treated flasks, dishes or multidishes and get a third case at no additional charge.



Qualifying Products for Promotion		
Description	Qty/Case	
Thermo Scientific Nunc EasyFlasks		
EasyFlask, Nunclon Delta-treated, Filter cap, 25cm <sup>2</sup> , 7mL working volume	10/200	
EasyFlask, Nunclon Delta-treated, Filter cap, 75cm <sup>2</sup> , 25mL working volume	5/100	
EasyFlask, Nunclon Delta-treated, Filter cap, 175cm <sup>2</sup> , 55mL working volume	5/30	
EasyFlask, Nunclon Delta-treated, Filter cap, 225cm <sup>2</sup> , 70mL working volume	5/30	
EasyFlask, Nunclon Delta-treated, Vent/close cap, 25cm <sup>2</sup> , 7mL working volume	10/200	
EasyFlask, Nunclon Delta-treated, Vent/close cap, 75cm <sup>2</sup> , 25mL working volume	5/100	
EasyFlask, Nunclon Delta-treated, Vent/close cap, 175cm <sup>2</sup> , 55mL working volume	5/30	
EasyFlask, Nunclon Delta-treated, Vent/close cap, 225cm <sup>2</sup> , 70mL working volume	5/30	
Thermo Scientific Nunc TripleFlasks		
TripleFlask, Nunclon Delta-treated, Vent/close cap, 500cm <sup>2</sup> , 200mL working volume	2/32	
TripleFlask, Nunclon Delta-treated, Filter cap, 500cm <sup>2</sup> , 200mL working volume	2/32	
TripleFlask with Bar Code, Nunclon Delta-treated, Filter cap, 500cm <sup>2</sup> , 200mL working volume	2/32	
Thermo Scientific Nunc Dishes		
35mm Dish, Nunclon Delta-treated, 8.8cm <sup>2</sup> , with lid; without airvent	10/500	
35mm Dish, Nunclon Delta-treated, 8.8cm <sup>2</sup> , with lid; with airvent	10/500	
35mm Dish, Nunclon Delta-treated, 8.8cm <sup>2</sup> , with lid; with airvent	20/500	
60mm Dish, Nunclon Delta-treated, 21.5cm <sup>2</sup> , with lid; without airvent	10/400	
60mm Dish, Nunclon Delta-treated, 21.5cm <sup>2</sup> , with lid; with airvent	10/400	
	Description   Thermo Scientific Nunc EasyFlasks   EasyFlask, Nunclon Delta-treated, Filter cap, 25cm², 7mL working volume   EasyFlask, Nunclon Delta-treated, Filter cap, 75cm², 25mL working volume   EasyFlask, Nunclon Delta-treated, Filter cap, 175cm², 55mL working volume   EasyFlask, Nunclon Delta-treated, Filter cap, 225cm², 70mL working volume   EasyFlask, Nunclon Delta-treated, Filter cap, 225cm², 70mL working volume   EasyFlask, Nunclon Delta-treated, Vent/close cap, 25cm², 7mL working volume   EasyFlask, Nunclon Delta-treated, Vent/close cap, 75cm², 25mL working volume   EasyFlask, Nunclon Delta-treated, Vent/close cap, 75cm², 25mL working volume   EasyFlask, Nunclon Delta-treated, Vent/close cap, 225cm², 70mL working volume   EasyFlask, Nunclon Delta-treated, Vent/close cap, 225cm², 70mL working volume   Thermo Scientific Nunc TripleFlasks   TripleFlask, Nunclon Delta-treated, Vent/close cap, 500cm², 200mL working volume   TripleFlask, Nunclon Delta-treated, Filter cap, 500cm², 200mL working volume   TripleFlask with Bar Code, Nunclon Delta-treated, Filter cap, 500cm², 200mL working volume   Thermo Scientific Nunc Dishes   35mm Dish, Nunclon Delta-treated, 8.8cm², with lid; without airvent   35mm Dish, Nunclon Delta-treated, 8.8cm², with lid; with airvent   35mm Dish, Nunclon Delta-treated, 8.8cm², with lid; with airvent   60mm Dish, Nunclon Delt	

(continued on next page)

Learn more about our flasks, dishes and multidishes

# engineered for superior and consistent cell growth

## **Thermo Scientific<sup>™</sup>Nunc<sup>™</sup>** Flasks, Dishes & Multidishes

Buy 2 cases of Nunclon™ Delta treated flasks, dishes or multidishes and get a third case at no additional charge.



Qualifying Products for Promotion			
Cat. No.	Description	Qty/Case	
172931	100mm Dish, Nunclon Delta-treated, 56.7cm <sup>2</sup> , with lid; with airvent	10/240	
150350	100mm Dish, Nunclon Delta-treated, 56.7 cm <sup>2</sup> , with lid; with airvent	150/150	
168381	150mm Dish, Nunclon Delta-treated, 145cm <sup>2</sup> , with lid; with airvent	10/80	
157150	150mm Dish, Nunclon Delta-treated, 150cm <sup>2</sup> , with lid; with airvent	10/120	
166508	245mm Dish, Nunclon Delta-treated, 500cm <sup>2</sup> , with lid; without airvent	4/16	
	Thermo Scientific Nunc Multidishes		
176740	4-well Multidish, Nunclon Delta-treated, 1.9cm <sup>2</sup> , round well; with lid	4/120	
167063	4-well Multidish, Nunclon Delta-treated, 21.8cm <sup>2</sup> , round well; with lid	10/100	
140675	6-well Multidish, Nunclon Delta-treated, 9.6cm <sup>2</sup> , round well; with lid	1/75	
140685	6-well Multidish, Nunclon Delta-treated, 9.6cm <sup>2</sup> , round well; with lid	5/85	
167064	8-well Multidish, Nunclon Delta-treated, 10.5cm <sup>2</sup> , rectangular well; with lid; sterile	10/100	
150628	12-well Multidish, Nunclon Delta-treated, 3.5cm <sup>2</sup> , round well; with lid; sterile	1/75	
142475	24-well Multidish, Nunclon Delta-treated, 1.8cm <sup>2</sup> , round well; with lid; sterile	1/75	
142485	24-well Multidish, Nunclon Delta-treated, 1.8cm <sup>2</sup> , round well; with lid; sterile	5/85	
150687	48-well Multidish, Nunclon Delta-treated, 1.1cm <sup>2</sup> , round well; with lid; sterile	1/75	
152640	48-well Multidish, Nunclon Delta-treated, 1.1cm <sup>2</sup> , round well; with lid; sterile	5/85	

Learn more about our flasks, dishes and multidishes

# superior performance in your cell-based assays

## **Thermo Scientific Nunc Edge and Cell Culture Plates**

Buy 1 case of cell culture plates and receive a second case at no additional charge.



Qualifying Products for Promotion		
Cat. No.	Description	Qty/Case
	Thermo Scientific Nunc Edge 96-well Plate	
167314	96 Well Plate PS, Edge Plate PS, Nunclon Delta Treated, with Lid, Clear F96	1/50
	96 Well Cell Culture Plates	
137101	96 Well Plate PS, Cell Culture Treated with Lid, Black F96 400µL Well	1/50
167008	96 Well Plate PS, Cell Culture Treated with Lid, Clear F96 400µL Well	1/50
163320	96 Well Plate PS, Cell Culture Treated with Lid, Clear U96 300µL Well	1/50
136101	96 Well Plate PS, Cell Culture Treated with Lid, White F96 400µL Well	1/50
168055	96 Well Plate PS, Cell Culture Treated without Lid, Clear F96 400µL Well	1/50
143761	96 Well Plate PS, Cell Culture Treated without Lid, Clear U96 300µL Well	1/50
165305	96 Well Plate PS, Opitcal Bottom Polymer Base, Cell Culture Treated with Lid, Black 400µL Well	10/30
165306	96 Well Plate PS, Opitcal Bottom Polymer Base, Cell Culture Treated with Lid, White 400 $\mu$ L Well	10/30
152037	96 Well Plate PS, Poly-D-Lysine Treated with Lid, Black/Clear Bottom 400µL Well	5/20
152039	96 Well Plate PS, Poly-D-Lysine Treated with Lid, Clear 300µL Well	5/20
152028	96 Well Plate PS, Poly-D-Lysine Treated with Lid White/Clear Bottom 400µL Well	5/20
165218	OmniTray PS, Nunclon Delta Treated with Lid, Clear 86x128mm	10/60

Learn more about our microplates for cell culture.

## your samples

### **Thermo Scientific Nalgene Cryogenic Tubes**



Buy any case of Nalgene cryogenic tubes and receive a Thermo Scientific<sup>™</sup> Nalgene<sup>™</sup> Mr. Frosty<sup>™</sup> of your choice at no additional charge.

#### **Qualifying Products for Promotion** Description Qty/Case Cat. No. **Nalgene Tubes** 25/500 5000-0012 Cryogenic Tube, 1.2mL, Sterile PP 5000-0020 Cryogenic Tube, 2.0mL, Sterile PP 25/500 5000-0050 Cryogenic Tube, 5.0mL, Sterile PP 10/250 5000-1012 System 100 Cryogenic Tube, 1.0mL, Sterile PP 25/5005000-1020 System 100 Cryogenic Tube, 1.5mL, Sterile PP 25/500 5001-0012 Cryogenic Tube with Linear Barcode, 1.2mL, Sterile PP 25/500 5001-0020 Cryogenic Tube with Linear Barcode, 2.0mL, Sterile PP 25/500 5001-0050 Cryogenic Tube with Linear Barcode, 5.0mL, Sterile PP 10/250 5011-1020 System 100 Cryogenic Tube with Barcode, 1.5mL, Sterile PP 25/500 5011-0012 Cryogenic Tube, Non-Sterile, Non-Printed, 1.2mL, Bulk pack PP 1000 5011-0020 Cryogenic Tube, Non-Sterile, Non-Printed, 2.0mL, Bulk pack PP 1000 5012-0012 Cryogenic Tube, Sterile, 1.2mL, Bulk pack PP 1000 5012-0020 Cryogenic Tube, Sterile, 2.0mL, Bulk pack PP 1000 Mr. Frosty 557100-200 1.0-2.0 Mr. Frosty 1 degree C Freezing Container, PC 1 18 89176-714 3.6 Mr. Frosty 1 degree C Freezing Container, PC 1 12 89176-716 4.5-5.0 Mr. Frosty 1 degree C Freezing Container, PC 1 12

Discover more crypreservation solutions for your cells.

#### flexible and efficient

sample storage and tracking



## **Thermo Scientific Universal Tubes**

Buy any case of Nunc Universal Tubes and receive a second case at no additional charge.

Qualifying Products for Promotion		
Cat. No.	Description	Qty/Case
374510	Nunc 2.0mL Tube, Internal Thread Universal Tube in Universal Latch Rack, 2D Barcoding, White Patch and Graduations	48/480
374511	Nunc 2.0mL Tube, Internal Thread Universal Tube in Universal Latch Rack, No Barcoding, White Patch and Graduations	48/480
374512	Nunc 2.0mL Tube, Internal Thread Universal Tube Bulk, 2D Barcoding, White Patch and Graduations	48/480
374513	Nunc 2.0mL Tube, Internal Thread Universal Tube Bulk, No Barcoding, White Patch and Graduations	48/480
374500	Nunc 1.8mL Tube, External Thread Universal Tube in Universal Latch Rack, 2D Barcoding, White Patch and Graduations	48/480
374501	Nunc 1.8mL Tube, External Thread Universal Tube Latch Rack, No Barcoding, White Patch and Graduations	48/480
374502	Nunc 1.8mL Tube, External Thread Universal Tube Bulk, 2D Barcoding, White Patch and Graduations	48/480
374503	Nunc 1.8mL Tube, External Thread Universal Tube Bulk, No Barcoding, White Patch and Graduations	48/480
374320	Nunc 5.0mL Low Profile Tube, External Thread Universal Tube Latch Rack, 2D Barcoding, White Patch and Graduations	24/240
374321	Nunc 5.0mL Low Profile Tube, External Thread Universal Tube Latch Rack, No Barcoding, White Patch and Graduations	24/240
374322	Nunc 5.0mL Low Profile Tube, External Thread Universal Tube Bulk, 2D Barcoding, White Patch and Graduations	24/240
374323	Nunc 5.0mL Low Profile Tube, External Thread Universal Tube Bulk, No Barcoding, White Patch and Graduations	24/240

Learn more about your Universal Tube options.

# ensuring reliable and repeatable microbiological quality control testing

#### Thermo Scientific<sup>™</sup> **Nalgene<sup>™</sup> QC Analytical Filters**

Buy 2 cases of QC Analytical Filters and receive the third case at no additional cost.



Qualifying Products for Promotion		
Cat. No.	Description	Qty/Case
130-4020	Filter Unit, 150 mL, 0.2µm CN Membrane, 47mm Diameter, Membrane Color/Grid: White/None	12/72
130-4045	Filter Unit, 150 mL, 0.45µm CN Membrane, 47mm Diameter, Membrane Color/Grid: White/Black	12/72
140-4045	Filter Funnel, 150 mL, 0.45µm CN Membrane, 47mm Diameter, Membrane Color/Grid: White Black	12/72
145-0020	Test Filter Funnel, 100 mL, 0.2µm CN Membrane, 47mm Diameter, Membrane Color/Grid: White/None	50
145-0045	Test Filter Funnel, 100 mL, 0.45µm CN Membrane, 47mm Diameter, Membrane Color/Grid: White/Black	50
145-2020	Test Filter Funnel, 250 mL, 0.2µm CN Membrane, 47mm Diameter, Membrane Color/Grid: White/None	50
145-2045	Test Filter Funnel, 250 mL, 0.45µm CN Membrane, 47mm Diameter, Membrane Color/Grid: White/Black	50
147-0045	Test Filter Funnel, 100 mL, 0.45µm CN Membrane, 47mm Diameter, Membrane Color/Grid: Grey/Black	50
147-2045	Test Filter Funnel, 250 mL, 0.45µm CN Membrane, 47mm Diameter, Membrane Color/Grid: Grey/Black	50

Find more information on all of our QC Analytical Filters.

Valid in the United States and Europe while supplies last for orders placed between February 2016 and June 30, 2016 and only in connection products identified above. Pricing in U.S. Dollars. Valid in U.S. only. All pricing information and values based on current manufacturer's suggested list price. "Savings" or "Discount" indicated herein represent the maximum difference between the standard manufacturer's suggested list price of the product and the manufacturer's suggested list price for the product under this promotion. Actual savings may vary, subject to dealer participation. Not valid in connection with any other offer or discount. Void where prohibited by law or company policy. Please note that items represented at no cost represent "in-kind discounts" and should be reported accordingly. In addition, the value of any specific offer(s) the customer may receive from the manufacturer under this program is a "discount or other reduction in price" to the customer under section 1128B(b)(3)(A) of the Social Security Act [42 U.S.C. 1320a-7b (b)(3)(a)]. Accordingly, the customer shall disclose this and any other discounts or other reduction in price received under this program under any state or federal program which provides cost or charge-based reimbursements to the customer for the products or services purchased under this program. Prices and specifications or products subject to change at any time. Not to be combined with any other offers. Please allow up to 8 weeks for delivery of promotional items. Please consult your local sales representative for details. Promotional offers begin February 1, 2016 and end June 30, 2016.

This promotion is not available to our current or potential healthcare customers. Healthcare customers include any healthcare practitioner, other healthcare provider or any individual or organization authorized to prescribe, dispense, purchase or influence the use of medical devices or supplies for clinical use. Healthcare customers include, without limitation, any "Health Care Practitioner" under the Massachusetts Marketing Code of Conduct Law (defined at 105 C.M.R §907.004) or a "Health Care Provider" under the Vermont Prescribed Products Gift Ban and Disclosure Law (defined at Vt. Stat. Ann. tit. 18, §4631a(8)). By participating in this promotion, you warrant that you are not a current or potential healthcare customer. Thermo Fisher reserves the right to determine your status in the event of a dispute. Customer acknowledges that this offer may include a discount or other price reduction that must be properly and accurately accounted for and reported by customer in accordance with all federal and state laws, including without limitation the federal anti-kickback law (42 C.F.R. § 1320a-7b(b)(3)(A)) and regulations thereunder (42 C.F.R. §1001.952(h)).

#### www.thermoscientific.com/cellculturepromos



**Thermo S C I E N T I F I C** A Thermo Fisher Scientific Brand

BRCCPROMOS 0116